AIRONE200[®]

LANÇAMENTOS!

**CÁLCIO ARSENAZO
MICROPROTE pirogalol**

REVISÃO: 05 (02/2009)

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TEMPERATURA: 37 °C

TEST NAME	ALB
DECIMAL	1
MEASURE UNIT	g/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	400
LINEARITY LIMIT	8.0
CONTAMINATING	YES
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	620
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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ALT/TGP UV

TEMPERATURA: 37 °C

TEST NAME	GPT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	800
MAXIMUM	2000
LINEARITY LIMIT	350
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 1746

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02

TEMPERATURA: 37 °C

TEST NAME	AMILAS
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0.0
MAXIMUM	500
LINEARITY LIMIT	2000
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	30
INCUBATION 2	0
LAG PHASE	15
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	3953

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AST/TGO UV

TEMPERATURA: 37 °C

TEST NAME	GOT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	800
MAXIMUM	2000
LINEARITY LIMIT	350
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 1746

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REVISÃO: 05 (02/2009)

TEST	CARD
Type	End Point.
Unit	mg/dL
Number of Decimals	1
Diluent	99 Diluen
Pre-washing	No
First Reagent Vol. (µL)	300
Second Reagent Vol. (µL)	0
Sample Volume (µL)	3.0
Lower Normal Value	8.8
Higer Normal Value	11.0
Systematic Dilution	1
Dilution for Urine	1
Low Reagent Blank	0.0000
High Reagent Blank	3.5000
Maximum Deviation	0.0100
Primary Wavelength	600 nm
High Wavelength	
Low Wavelength	
Systematic Serum Blank	No
T Tolerance Factor	0.0150
Antigen Excess Test (%)	No
Slope	1.00
Intercept	0.0
Reaction Delay	02:00

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CK - NAC

TEMPERATURA: 37°C

TEST NAME	CK
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	600
LINEARITY LIMIT	1700
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	120
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	8095

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REVISÃO: 05 (02/2009)

TEMPERATURA: 37 °C

TEST NAME	CLOC
DECIMAL	1
MEASURE UNIT	mmol/L
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	130
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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COLESTEROL 250

TEMPERATURA: 37 °C

TEST NAME	COLT
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	500
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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REVISÃO: 05 (02/2009)

TEMPERATURA: 37°C

TEST NAME	COLT
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	800
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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COLESTEROL HDL

TEMPERATURA: 37°C

TEST NAME	CHDL
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	10
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	400
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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REVISÃO: 05 (02/2009)

TEMPERATURA: 37 °C

TEST NAME	CREAT
DECIMAL	2
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	500
LINEARITY LIMIT	10.0
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	20
INCUBATION 2	0
LAG PHASE	10
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0

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DHL - UV

TEMPERATURA: 37 °C

TEST NAME	LDH
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	800
MAXIMUM	2000
LINEARITY LIMIT	1257
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	30
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	- 8109

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TEMPERATURA: 37 °C

TEST NAME	FALC
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	2000
LINEARITY LIMIT	690
CONTAMINATING	YES
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	60
INCUBATION 2	0
LAG PHASE	15
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	2764

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FOSFATO - UV

TEMPERATURA: 37 °C

TEST NAME	FOSF
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1500
LINEARITY LIMIT	20
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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REVISÃO: 05 (02/2009)

TEMPERATURA: 37 °C

TEST NAME	GGT
DECIMAL	0
MEASURE UNIT	U/L
REACTION TYPE	KINETIC
SAMPLE VOLUME (µL)	40
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	1000
LINEARITY LIMIT	300
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	405
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	30
INCUBATION 2	0
LAG PHASE	10
MEASURE	30
MEASUREMENT TYPE	
CALIBRATE	FACTOR
FACTOR	1111

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GLICOSE ENZ. LÍQUIDA

TEMPERATURA: 37 °C

TEST NAME	GLIC
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	600
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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REVISÃO: 05 (02/2009)

TEMPERATURA: 37°C

TEST NAME	PROT.TOT
DECIMAL	1
MEASURE UNIT	g/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	250
LINEARITY LIMIT	12.0
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	546
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	300
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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TRIGLICÉRIDES 120

TEMPERATURA: 37°C

TEST NAME	TGR
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	1000
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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REVISÃO: 05 (02/2009)

TEMPERATURA: 37 °C

TEST NAME	TGR
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	800
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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URATO 160

TEMPERATURA: 37 °C

TEST NAME	AU
DECIMAL	1
MEASURE UNIT	mg/dL
REACTION TYPE	END POINT
SAMPLE VOLUME (µL)	8
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	0
MAXIMUM	300
LINEARITY LIMIT	14.0
CONTAMINATING	NO
REAGENT BLANKING	YES
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	510
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	600
INCUBATION 2	0
LAG PHASE	3
MEASURE	1
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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URÉIA - UV

TEMPERATURA: 37°C

TEST NAME	UREUV
DECIMAL	0
MEASURE UNIT	mg/dL
REACTION TYPE	FIXED TIME
SAMPLE VOLUME (µL)	4
REAGENT VOLUME (µL)	
REAGENT 1	400
REAGENT 2	0
ABSORBANCE RANGE (mABS)	
MINIMUM	1000
MAXIMUM	2000
LINEARITY LIMIT	200
CONTAMINATING	NO
REAGENT BLANKING	NO
DIFFERENTIAL	NO
FILTERS (nm)	
FILTER 1	340
FILTER 2	NONE
TIMES (SEC)	
MIX 1	0.00
MIX 2	0.00
INCUBATION 1	20
INCUBATION 2	0
LAG PHASE	10
MEASURE	60
MEASUREMENT TYPE	
CALIBRATE	STANDARD
FACTOR	0.0

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