AUTOHUMALYSER 900 S[®]

LANÇAMENTOS!

**CÁLCIO ARSENAZO
MICROPROTE pirogalol**

REVISÃO: 04 (02/2009)

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TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | ALBUMINA |
| DECIMAL | 0 |
| MEASURE UNIT | g/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 400 |
| LINEARITY LIMIT | 8.0 |
| CONTAMINATING | YES |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 620 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 120 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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ALT/TGP UV

TEMPERATURA: 37°C

| | |
|--------------------------------|---------|
| TESTNAME | ALT/GPT |
| DECIMAL | 0 |
| MEASURE UNIT | U/L |
| REACTION TYPE | KINETIC |
| SAMPLE VOLUME (µL) | 40 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 1000 |
| MAXIMUM | 2000 |
| LINEARITY LIMIT | 350 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 340 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 60 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 30 |
| MEASUREMENT TYPE | |
| CALIBRATE | FACTOR |
| FACTOR | -1746 |

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REVISÃO: 04 (02/2009)

TEMPERATURA 37°C

| | |
|--------------------------------|---------|
| TESTNAME | AMICNPG |
| DECIMAL | 0 |
| MEASURE UNIT | U/l |
| REACTION TYPE | KINETIC |
| SAMPLE VOLUME (µL) | 8 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0.0 |
| MAXIMUM | 500 |
| LINEARITY LIMIT | 2000 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 405 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 60 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 30 |
| MEASUREMENT TYPE | |
| CALIBRATE | FACTOR |
| FACTOR | 3953 |

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AST/TGO UV

TEMPERATURA: 37°C

| | |
|--------------------------------|---------|
| TESTNAME | AST/GOT |
| DECIMAL | 0 |
| MEASURE UNIT | U/L |
| REACTION TYPE | KINETIC |
| SAMPLE VOLUME (µL) | 40 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 1000 |
| MAXIMUM | 2000 |
| LINEARITY LIMIT | 350 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 340 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 60 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 30 |
| MEASUREMENT TYPE | |
| CALIBRATE | FACTOR |
| FACTOR | -1746 |

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REVISÃO: 04 (02/2009)

TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | CÁLCIO |
| DECIMAL | 1 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 3 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 300 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 1000 |
| LINEARITY LIMIT | 20 |
| CONTAMINATING | YES |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 620 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 120 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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CK - NAC

TEMPERATURA: 37°C

| | |
|--------------------------------|---------|
| TESTNAME | CK |
| DECIMAL | 0 |
| MEASURE UNIT | U/L |
| REACTION TYPE | KINETIC |
| SAMPLE VOLUME (µL) | 8 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 500 |
| LINEARITY LIMIT | 1700 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 340 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.40 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 120 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 30 |
| MEASUREMENT TYPE | |
| CALIBRATE | FACTOR |
| FACTOR | 8095 |

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TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | CLORETOS |
| DECIMAL | 0 |
| MEASURE UNIT | mEq/L |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | |
| LINEARITY LIMIT | 125 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 510 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0,30 |
| MIX 2 | |
| INCUBATION 1 | 120 |
| INCUBATION 2 | 0 |
| LAG PHASE | 2 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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COLESTEROL 250

TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | COLT |
| DECIMAL | 0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 300 |
| LINEARITY LIMIT | 500 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 510 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.40 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | COLT |
| DECIMAL | 0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 300 |
| LINEARITY LIMIT | 800 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 510 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0,00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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COLESTEROL HDL

TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | CHDL |
| DECIMAL | 0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 10 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 300 |
| LINEARITY LIMIT | 400 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 510 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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REVISÃO: 04 (02/2009)

TEMPERATURA: 37°C

| | |
|--------------------------------|------------|
| TESTNAME | CREA |
| DECIMAL | 0.0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | FIXED TIME |
| SAMPLE VOLUME (µL) | 40 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 500 |
| LINEARITY LIMIT | 10.0 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | |
| FILTERS (nm) | |
| FILTER 1 | 510 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 30 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 60 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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DHL - UV

TEMPERATURA: 37°C

| | |
|--------------------------------|---------|
| TEST NAME | LDH |
| DECIMAL | 0 |
| MEASURE UNIT | U/L |
| REACTION TYPE | KINETIC |
| SAMPLE VOLUME (µL) | 8 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 800 |
| MAXIMUM | 2000 |
| LINEARITY LIMIT | 1257 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 340 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 60 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 30 |
| MEASUREMENT TYPE | |
| CALIBRATE | FACTOR |
| FACTOR | - 8109 |

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TEMPERATURA: 37°C

| | |
|--------------------------------|---------|
| TESTNAME | FALC |
| DECIMAL | 0 |
| MEASURE UNIT | U/L |
| REACTION TYPE | KINETIC |
| SAMPLE VOLUME (µL) | 8 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 2000 |
| LINEARITY LIMIT | 690 |
| CONTAMINATING | YES |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 405 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 60 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 30 |
| MEASUREMENT TYPE | |
| CALIBRATE | FACTOR |
| FACTOR | 2764 |

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FOSFATO - UV

TEMPERATURA: 37°C

| | |
|--------------------------------|----------|
| TESTNAME | FÓSFORO |
| DECIMAL | 0.0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END |
| POINT | |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 1500 |
| LINEARITY LIMIT | 20 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 340 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 300 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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TEMPERATURA: 37°C

| | |
|--------------------------------|---------|
| TESTNAME | GGT |
| DECIMAL | 0 |
| MEASURE UNIT | U/L |
| REACTION TYPE | KINETIC |
| SAMPLE VOLUME (µL) | 40 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 750 |
| LINEARITY LIMIT | 300 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 405 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 60 |
| INCUBATION 2 | 0 |
| LAG PHASE | 15 |
| MEASURE | 30 |
| MEASUREMENT TYPE | |
| CALIBRATE | FACTOR |
| FACTOR | 1111 |

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GLICOSE ENZ. LÍQUIDA

TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | GLIC |
| DECIMAL | 0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 300 |
| LINEARITY LIMIT | 600 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 510 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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MICROPROTE pirogalol

TEMPERATURA: 37°C

| | |
|--------------------------------|-------------|
| TESTNAME | PROTEINÚRIA |
| DECIMAL | 1 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 6 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 300 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 400 |
| LINEARITY LIMIT | 300 |
| CONTAMINATING | YES |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 620 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TEST NAME | PROT.TOT |
| DECIMAL | 1 |
| MEASURE UNIT | g/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 8 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 250 |
| LINEARITY LIMIT | 12.0 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 546 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 300 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0.0 |

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TRIGLICÉRIDES 120

TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | TGR |
| DECIMAL | 0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 300 |
| LINEARITY LIMIT | 1000 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 546 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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REVISÃO: 04 (02/2009)

TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | TGR |
| DECIMAL | 0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 300 |
| LINEARITY LIMIT | 800 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 546 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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URATO 160

TEMPERATURA: 37°C

| | |
|--------------------------------|-----------|
| TESTNAME | AU |
| DECIMAL | 0.0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | END POINT |
| SAMPLE VOLUME (µL) | 8 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 0 |
| MAXIMUM | 300 |
| LINEARITY LIMIT | 14.0 |
| CONTAMINATING | NO |
| REAGENT BLANKING | YES |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 510 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 600 |
| INCUBATION 2 | 0 |
| LAG PHASE | 3 |
| MEASURE | 1 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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URÉIA - UV

TEMPERATURA: 37°C

| | |
|--------------------------------|------------|
| TEST NAME | UREUV |
| DECIMAL | 0 |
| MEASURE UNIT | mg/dL |
| REACTION TYPE | FIXED TIME |
| SAMPLE VOLUME (µL) | 4 |
| REAGENT VOLUME (µL) | |
| REAGENT 1 | 400 |
| REAGENT 2 | 0 |
| ABSORBANCE RANGE (mABS) | |
| MINIMUM | 1000 |
| MAXIMUM | 2000 |
| LINEARITY LIMIT | 200 |
| CONTAMINATING | NO |
| REAGENT BLANKING | NO |
| DIFFERENTIAL | NO |
| FILTERS (nm) | |
| FILTER 1 | 340 |
| FILTER 2 | NONE |
| TIMES (SEC) | |
| MIX 1 | 0.00 |
| MIX 2 | 0.00 |
| INCUBATION 1 | 30 |
| INCUBATION 2 | 0 |
| LAG PHASE | 10 |
| MEASURE | 60 |
| MEASUREMENT TYPE | |
| CALIBRATE | STANDARD |
| FACTOR | 0 |

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